

## Unistat® 425

### Cooling a DDPS 2-litre jacketed glass reactor to $T_{min}$

#### Requirement

This case study is to find out the minimum temperature that a Unistat 425 can cool the jacket of a 2-litre jacketed glass reactor and the resultant process temperature.

#### Method

The Unistat 425 is connected to the 2-litre DDPS glass reactor using two insulated metal 1-metre hoses. The reactor is filled with 1.5 litre of "M90.055.03", a silicon based HTF

#### Results

The jacket is cooled to  $-40\text{ }^{\circ}\text{C}$  in around 33 minutes. The process temperature ramp rate slows as the  $\Delta T$  narrows and has reached  $-39\text{ }^{\circ}\text{C}$  when the test is stopped.

#### Setup details

Unistat® 425 & DDPS reactor

Temperature range:  $-40\text{...}250\text{ }^{\circ}\text{C}$   
 Cooling power:  $2.5\text{ kW @ }0\text{ }^{\circ}\text{C}$   
 $1.8\text{ kW @ }-20\text{ }^{\circ}\text{C}$

Heating power:  $2.0\text{ kW}$   
 Hoses:  $2 \times 1\text{ m; M24} \times 1.5$   
 (#9325)

HTF: DW-Therm (#6479)  
 Reactor: 2-litre jacketed glass reactor

Reactor contents:  $1.5\text{ litre M90.055.03}$   
 (#6259)

Reactor stirrer speed:  $150\text{ rpm}$   
 Control: internal

