

Ministat® 240

Ministat® 240 cycling a 20-liter glass jacketed reactor

Requirement

This Case Study demonstrates the minimum achievable process temperature and the temperature control capabilities when Ministat 240 is connected with a Buchi Glas Üster 20-liter reactor.

Method

The 20-liter Buchi Glas Üster reactor was connected to Ministat® 240 using 1-meter metal insulated hoses. The thermofluid used in the system was "DW-Therm". "Process" control was carried out via a Pt100 sensor located in the process mass. Stirrer speed was set to 100 rpm.

Setup details

- Temperature range: -45°C...+200°C
- Cooling power: 0.60 kW @ +20°C
- 0.55 kW @ 0°C
- 0.35 kW @ -20°C
- Heating power: 2.0 kW
- Hoses: 2*1 m metal insulated
- HTF: DW-Therm
- Reactor: Buchi Glas Üster 20-liter glass jacketed
- Reactor content: 15 l DW-Therm
- Stirrer speed: 100 rpm
- Control: process
- Amb. temperature: +24°C

Video clip


<https://youtu.be/7nlvU3VKqBQ>



Results

1. Lowest achievable temperature (Tmin):

As the graphic shows that the minimum achievable Process Temperature was -17.6°C.



2. Performance:

The graphic shows the speed, accuracy and stability as each new set point is reached.

Start T	End T	Approximate time	Av. Ramp Rate
+20°C	0°C	69 minutes	0.3 K/min
0°C	+100°C	95 minutes	1.1 K/min
+100°C	+20°C	134 minutes	0.6 K/min

