

# Unistat® 705w

**Cooling a Buchi Glas Uster 3-litre metal reactor from 180 °C to 20 °C**

**Requirement**

This case study shows the performance of a Unistat 705w cooling a Buchi Glas Uster 3-litre un-insulated metal pressure reactor from 180 °C to 20 °C under “process” control.

**Method**

The Unistat and reactor are connected using two 1-metre insulated metal hoses. The reactor is filled with 2.25 litre of “M90.055.03”, a Huber supplied silicon based HTF.

**Results**

It can be seen from the narrow  $\Delta T$  between “internal” (jacket) and the process temperatures that thermal transfer in this reactor is very good. Total time to ramp the process through 160 K (180 °C to 20 °C) is approximately 65 minutes.

**Setup details**

Unistat® 705w & Buchi Glas Uster reactor

- Temperature range: -75...250 °C
- Cooling Power: 0.6 kW @ 250...100 °C  
0.65 kW @ 0 °C  
0.6 kW @ -20...-40 °C  
0.3 kW @ -60 °C
- Heating Power: 1.5 kW / 3 kW
- Pump speed: 3500 rpm
- Hoses: 2x1 m; M24x1.5 (#9325)
- HTF: DW-Therm (#6479)
- Reactor: 3-litre un-insulated metal pressure reactor
- Reactor contents: 2.25 litre M90.055.03 (#6259)
- Reactor stirrer speed: 200 rpm
- Control: process

