

External pressure pump. Bearing out of the liquid circuit. All parts in the liquid circuit and housing made of stainless steel. Operation either manual (on/off) or external via 3-line-connector. Cable (Option) for the potential free contact of Unistat Tango and the big power Unistats, all Compatible Control Thermostats and Unichillers.

Technical data according to DIN 12876

Operating temperature range	-120...300 °C	
max. delivery of pressure pump	70 l/min	
max. delivery pressure	1.5 bar	Order-No.: 527.0021
Pump connection	M38x1,5 male	
Overall dimensions WxDxH **	457x201x320 mm	
Net weight	20 kg	
Power supply (3 Phase)	400V 3-N 50Hz	
max. current (3 Phase)	2 A	
Fuse (3 phase)	3x4 A	
min. ambient temperature	5 °C	
max. ambient temperature	32 °C	
from Serial-No.:	193694	1.0/14

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions. Illustrations can deviate from the original.

Output data valid for: Room temperature 20° C

Standard delivery conditions - Power cable configuration:

1. Single-phase devices (230V/115V) -> with cable and plug
2. Three-phase devices with current consumption less than 63A -> with cable, without plug
3. Three-phase devices with current consumption greater than 63A -> without cable, without plug

** Please respect space requirements. See operating conditions at www.huber-online.com